




Metallo-beta-lactamase type 2 Antibody

Product Code	CSB-PA346782ZA01SYN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52699
Immunogen	Recombinant <i>Serratia marcescens</i> Metallo-beta-lactamase type 2 protein (19-246aa)
Species Reactivity	<i>Serratia marcescens</i>
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	<i>Serratia marcescens</i>
Research Area	Others
Target Names	N/A

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Metallo-beta-lactamase type 2
antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 32.3 kDa
Observed band size: 33 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.