



Dct Antibody

Product Code	CSB-PA006562ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29812
Immunogen	Recombinant Mus musculus Dct protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-dopachrome tautomerase (DCT) (DT) (EC 5.3.3.12) (DOPAchrome conversion factor) (DOPAchrome isomerase) (DOPAchrome oxidoreductase) (L-dopachrome Delta-isomerase) (SLATY locus protein) (Tyrosinase-related protein 2) (TRP-2) (TRP2), Dct, Tyrp-2 Tyrp2
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Dct
Product Modify	Unmodified