



UL131 Antibody

Product Code	CSB-PA323723ZA01HWV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16773
Immunogen	Recombinant Human cytomegalovirus UL131 protein
Raised In	Rabbit
Species Reactivity	Human cytomegalovirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein UL131, UL131
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (Human herpesvirus 5)
Lead Time	Made-to-order (12-14 weeks)
Target Names	UL131
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.