



ORF2 Antibody

Product Code	CSB-PA897430ZA01HHAM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9YLQ9
Immunogen	Recombinant Hepatitis E virus genotype 3 ORF2 protein
Raised In	Rabbit
Species Reactivity	Hepatitis E virus genotype 3
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Capsid protein (Protein ORF2) (pORF2)
Immunogen Species	Hepatitis E virus genotype 3 (isolate Human/United States/US2) (HEV-3)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ORF2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.