



CD40LG Antibody

Product Code	CSB-PA004937ZA01MOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63304
Immunogen	Recombinant Macaca mulatta CD40LG protein
Raised In	Rabbit
Species Reactivity	Macaca mulatta
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD40 ligand (CD40-L) (Tumor necrosis factor ligand superfamily member 5) (CD antigen CD154) [Cleaved into: CD40 ligand, membrane form; CD40 ligand, soluble form], CD40LG, CD40L TNFSF5
Immunogen Species	Macaca mulatta (Rhesus macaque)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CD40LG
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.