



UL111A Antibody

Product Code	CSB-PA517245ZA01HWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F5HC71
Immunogen	Recombinant Human cytomegalovirus UL111A protein
Raised In	Rabbit
Species Reactivity	Human cytomegalovirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Viral interleukin-10 homolog (cmvIL-10) (vIL-10), UL111A
Immunogen Species	Human cytomegalovirus (strain Merlin) (HHV-5) (Human herpesvirus 5)
Lead Time	Made-to-order (12-14 weeks)
Target Names	UL111A
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.