



IL2 Antibody

Product Code	CSB-PA011629ZA01MOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68291
Immunogen	Recombinant Macaca mulatta IL2 protein
Raised In	Rabbit
Species Reactivity	Macaca mulatta
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-2 (IL-2) (T-cell growth factor) (TCGF), IL2
Immunogen Species	Macaca mulatta (Rhesus macaque)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IL2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.