



# FIB1 Antibody

<b>Product Code</b>	CSB-PA861762ZA01DOA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9FEF8
<b>Immunogen</b>	Recombinant Arabidopsis thaliana FIB1 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Probable mediator of RNA polymerase II transcription subunit 36b (Histone-glutamine methyltransferase) (SKP1-interacting partner 7) (rRNA 2'-O-methyltransferase fibrillar 1) (EC 2.1.1.-), MED36B, FBR1 FIB1 MED36_2 SKIP7
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	FIB1
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.