



mdh Antibody

Product Code	CSB-PA330043ZA01BAU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31005
Immunogen	Recombinant Bacillus methanolicus mdh protein
Raised In	Rabbit
Species Reactivity	Bacillus methanolicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-dependent methanol dehydrogenase (MDH) (MEDH) (EC 1.1.1.244) (Type 3 alcohol dehydrogenase), mdh
Immunogen Species	Bacillus methanolicus
Lead Time	Made-to-order (12-14 weeks)
Target Names	mdh
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.