



sod1 Antibody

Product Code	CSB-PA022397ZA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O73872
Immunogen	Recombinant Danio rerio sod1 protein
Raised In	Rabbit
Species Reactivity	Danio rerio
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Superoxide dismutase [Cu-Zn] (EC 1.15.1.1), sod1, cuzn
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Lead Time	Made-to-order (12-14 weeks)
Target Names	sod1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.