



rnj Antibody

Product Code	CSB-PA344938ZA01SSQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54123
Immunogen	Recombinant Synechocystis sp. rnj protein
Raised In	Rabbit
Species Reactivity	Synechocystis sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease J (RNase J) (EC 3.1.-.-), rnj
Immunogen Species	Synechocystis sp. (strain PCC 6803 / Kazusa)
Lead Time	Made-to-order (12-14 weeks)
Target Names	rnj
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.