



mppA Antibody

Product Code	CSB-PA304378ZA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P77348
Immunogen	Recombinant Escherichia coli mppA protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Periplasmic murein peptide-binding protein, mppA, ynaH
Immunogen Species	Escherichia coli (strain K12)
Lead Time	Made-to-order (12-14 weeks)
Target Names	mppA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.