



Igfbp1 Antibody

Product Code	CSB-PA011095ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21743
Immunogen	Recombinant Rat Insulin-like growth factor-binding protein 1 (Igfbp1) (26-272AA)
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin-like growth factor-binding protein 1 (IBP-1) (IGF-binding protein 1) (IGFBP-1), Igfbp1, Igfbp-1
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order
Target Names	Igfbp1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.