



Masp2 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA856307ZA01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q91WP0 |
| Immunogen | Recombinant Mus musculus Masp2 protein |
| Raised In | Rabbit |
| Species Reactivity | Mus musculus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Mannan-binding lectin serine protease 2 (EC 3.4.21.104) (MBL-associated serine protease 2) (Mannose-binding protein-associated serine protease 2) (MASP-2) [Cleaved into: Mannan-binding lectin serine protease 2 A chain; Mannan-binding lectin serine protease 2 B chain], Masp2 |
| Immunogen Species | Mus musculus (Mouse) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | Masp2 |
| Product Modify | Unmodified |