



dhaK Antibody

Product Code	CSB-PA338209ZA01DTO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45510
Immunogen	Recombinant Citrobacter freundii dhaK protein
Raised In	Rabbit
Species Reactivity	Citrobacter freundii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydroxyacetone kinase (DHA kinase) (EC 2.7.1.29) (Glycerone kinase), dhaK
Immunogen Species	Citrobacter freundii
Lead Time	Made-to-order (12-14 weeks)
Gene Names	dhaK
Product Modify	Unmodified