



ALS3 Antibody

Product Code	CSB-PA527942ZA01CZF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O74623
Immunogen	Recombinant Candida albicans ALS3 protein
Raised In	Rabbit
Species Reactivity	Candida albicans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Agglutinin-like protein 3 (3D9 antigen) (Adhesin 3), ALS3, ALD8
Immunogen Species	Candida albicans (Yeast)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ALS3
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.