



Putative serine protease inhibitor Antibody

Product Code	CSB-PA21189ZA01WBS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F6KME1
Immunogen	Recombinant Wuchereria bancrofti Putative serine protease inhibitor protein
Raised In	Rabbit
Species Reactivity	Wuchereria bancrofti
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative serine protease inhibitor
Immunogen Species	Wuchereria bancrofti
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified