



adk Antibody

Product Code	CSB-PA769467ZA01SZB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q83M40
Immunogen	Recombinant Shigella flexneri adk protein
Raised In	Rabbit
Species Reactivity	Shigella flexneri
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenylate kinase (AK) (EC 2.7.4.3) (ATP-AMP transphosphorylase) (ATP:AMP phosphotransferase) (Adenylate monophosphate kinase), adk
Immunogen Species	Shigella flexneri
Lead Time	Made-to-order (12-14 weeks)
Gene Names	adk
Product Modify	Unmodified