



amiB2 Antibody

Product Code	CSB-PA351872ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WQ97
Immunogen	Recombinant Mycobacterium tuberculosis amiB2 protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative amidase AmiB2 (EC 3.5.1.4), amiB2
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	amiB2
Product Modify	Unmodified