



fanC Antibody

Product Code	CSB-PA322992ZA01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18103
Immunogen	Recombinant Escherichia coli fanC protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	K99 fimbrial protein, fanC
Immunogen Species	Escherichia coli
Lead Time	Made-to-order (12-14 weeks)
Gene Names	fanC
Product Modify	Unmodified