



PTPRC Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA019049ZA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P08575 |
| Immunogen | Recombinant Homo sapiens PTPRC protein |
| Raised In | Rabbit |
| Species Reactivity | Homo sapiens |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Receptor-type tyrosine-protein phosphatase C (EC 3.1.3.48) (Leukocyte common antigen) (L-CA) (T200) (CD antigen CD45), PTPRC, CD45 |
| Immunogen Species | Homo sapiens (Human) |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | PTPRC |
| Product Modify | Unmodified |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |