



Dct Antibody

Product Code	CSB-PA006562ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29812
Immunogen	Recombinant Mouse L-dopachrome tautomerase (Dct), partial (24-472AA)
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-dopachrome tautomerase (DCT) (DT) (EC 5.3.3.12) (DOPAchrome conversion factor) (DOPAchrome isomerase) (DOPAchrome oxidoreductase) (L-dopachrome Delta-isomerase) (SLATY locus protein) (Tyrosinase-related protein 2) (TRP-2) (TRP2), Dct, Tyrp-2 Tyrp2
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order
Target Names	Dct
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.