



ADH2 Antibody

Product Code	CSB-PA842610ZA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96533
Immunogen	Recombinant Arabidopsis thaliana Alcohol dehydrogenase class-3 (ADH2) (2-379AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alcohol dehydrogenase class-3 (EC 1.1.1.1) (Alcohol dehydrogenase class-III) (Glutathione-dependent formaldehyde dehydrogenase) (FALDH) (FDH) (GSH-FDH) (EC 1.1.1.-) (S-(hydroxymethyl)glutathione dehydrogenase) (EC 1.1.1.284), ADH2, ADHIII FDH1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order
Target Names	ADH2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.