



# IpxC Antibody

<b>Product Code</b>	CSB-PA409457ZA01PZL
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	A6VB79
<b>Immunogen</b>	Recombinant Pseudomonas aeruginosa IpxC protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Pseudomonas aeruginosa
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	UDP-3-O-acyl-N-acetylglucosamine deacetylase (UDP-3-O-acyl-GlcNAc deacetylase) (EC 3.5.1.108) (UDP-3-O-[R-3-hydroxymyristoyl]-N-acetylglucosamine deacetylase), IpxC
<b>Immunogen Species</b>	Pseudomonas aeruginosa (strain PA7)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	IpxC
<b>Product Modify</b>	Unmodified