



LCE3A Antibody

Product Code	CSB-PA716360ZA01HUb0
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TA76
Immunogen	Recombinant Homo sapiens LCE3A protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Late cornified envelope protein 3A (Late envelope protein 13), LCE3A, LEP13
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	LCE3A
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.