



Hnrnpa2b1 Antibody

Product Code	CSB-PA010602ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O88569
Immunogen	Recombinant Mouse Heterogeneous nuclear ribonucleoproteins A2/B1 (Hnrnpa2b1) (1-353AA)
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A2/B1), Hnrnpa2b1, Hnrpa2b1
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order
Target Names	Hnrnpa2b1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.