



Hnrnpa2b1 Antibody

| | |
|---------------------|----------------------------------------------------------------------------------------|
| Product Code | CSB-PA010602ZA01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O88569 |
| Immunogen | Recombinant Mouse Heterogeneous nuclear ribonucleoproteins A2/B1 (Hnrnpa2b1) (1-353AA) |
| Raised In | Rabbit |
| Species Reactivity | Mus musculus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A2/B1), Hnrnpa2b1, Hnrpa2b1 |
| Immunogen Species | Mus musculus (Mouse) |
| Lead Time | Made-to-order |
| Gene Names | Hnrnpa2b1 |
| Product Modify | Unmodified |