



IaaA Antibody

Product Code	CSB-PA748477ZA01PHAЕ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76KX0
Immunogen	Recombinant Pseudomonas azotoformans IaaA protein
Raised In	Rabbit
Species Reactivity	Pseudomonas azotoformans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-amino acid amidase (EC 3.5.1.101), IaaA
Immunogen Species	Pseudomonas azotoformans
Lead Time	Made-to-order (12-14 weeks)
Target Names	IaaA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.