



subC Antibody

Product Code	CSB-PA365470ZA01BQT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00780
Immunogen	Recombinant Bacillus licheniformis subC protein
Raised In	Rabbit
Species Reactivity	Bacillus licheniformis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Subtilisin Carlsberg (EC 3.4.21.62), apr
Immunogen Species	Bacillus licheniformis
Lead Time	Made-to-order (12-14 weeks)
Target Names	subC
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.