



Mylk Antibody

Product Code	CSB-PA320904ZA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11799
Immunogen	Recombinant Gallus gallus Mylk protein
Raised In	Rabbit
Species Reactivity	Gallus gallus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin light chain kinase, smooth muscle (MLCK) (EC 2.7.11.18) (Telokin) [Cleaved into: Myosin light chain kinase, smooth muscle, deglutamylated form], Mylk
Immunogen Species	Gallus gallus (Chicken)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Mylk
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.