



2S albumin Antibody

Product Code	CSB-PA308526ZA01TLH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P86783
Immunogen	Recombinant Taraxacum officinale 2S albumin protein
Raised In	Rabbit
Species Reactivity	Taraxacum officinale
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	2S albumin (To-A1) [Cleaved into: 2S albumin small chain; 2S albumin large chain] (Fragments)
Immunogen Species	Taraxacum officinale (Common dandelion) (Leontodon taraxacum)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.