



Tpsab1 Antibody

Product Code	CSB-PA024127ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27435
Immunogen	Recombinant Rattus norvegicus Tpsab1 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptase (EC 3.4.21.59) (Mast cell protease 7) (rMCP-7) (Tryptase alpha/beta-1) (Tryptase, skin), Tpsab1, Mcp7 Mcpt7 Tpsb1
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Tpsab1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.