



fabB Antibody

Product Code	CSB-PA359247ZA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A953
Immunogen	Recombinant Escherichia coli fabB protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-oxoacyl-[acyl-carrier-protein] synthase 1 (EC 2.3.1.41) (3-oxoacyl-[acyl-carrier-protein] synthase I) (Beta-ketoacyl-ACP synthase I) (KAS I), fabB, fabC
Immunogen Species	Escherichia coli (strain K12)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	fabB
Product Modify	Unmodified