



# fabB Antibody

<b>Product Code</b>	CSB-PA359247ZA01ENV
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0A953
<b>Immunogen</b>	Recombinant Escherichia coli 3-oxoacyl-[acyl-carrier-protein] synthase 1 (fabB) (1-406AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Escherichia coli
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Protein A/G
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	3-oxoacyl-[acyl-carrier-protein] synthase 1 (EC 2.3.1.41) (3-oxoacyl-[acyl-carrier-protein] synthase I) (Beta-ketoacyl-ACP synthase I) (KAS I), fabB, fabC
<b>Species</b>	Escherichia coli (strain K12)
<b>Lead Time</b>	Made-to-order
<b>Gene Names</b>	fabB
<b>Product Modify</b>	Unmodified