



hfb2 Antibody

Product Code	CSB-PA304437ZA01HYE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79073
Immunogen	Recombinant Hypocrea jecorina Hydrophobin-2 (hfb2) (16-86AA)
Raised In	Rabbit
Species Reactivity	Hypocrea jecorina
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydrophobin-2 (Hydrophobin II) (HFBII), hfb2
Immunogen Species	Hypocrea jecorina (Trichoderma reesei)
Lead Time	Made-to-order
Gene Names	hfb2
Product Modify	Unmodified