



relG Antibody

Product Code	CSB-PA521896ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O33348
Immunogen	Recombinant Mycobacterium tuberculosis relG protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toxin RelG (EC 3.1.-.-) (Putative endoribonuclease RelG), relG, relE2
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Lead Time	Made-to-order (12-14 weeks)
Target Names	relG
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.