



IL5 Antibody

Product Code	CSB-PA863067ZA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9MYM5
Immunogen	Recombinant Sus scrofa IL5 protein
Raised In	Rabbit
Species Reactivity	Sus scrofa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-5 (IL-5) (Eosinophil differentiation factor) (T-cell replacing factor) (TRF), IL5
Immunogen Species	Sus scrofa (Pig)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IL5
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.