



tsaE Antibody

Product Code	CSB-PA360227ZA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AF67
Immunogen	A synthesized peptide derived from Escherichia coli tsaE protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, Dot blot (guaranteed specific to the immunizing peptide)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA threonylcarbamoyladenosine biosynthesis protein TsaE (t(6)A37 threonylcarbamoyladenosine biosynthesis protein TsaE), tsaE, yjeE
Immunogen Species	Escherichia coli (strain K12)
Lead Time	Made-to-order
Target Names	tsaE
Product Modify	①. 100ug * KLH-conjugated peptide paired with 100ug * BSA-conjugated peptide (positive control); ②. 1ml*pre-immune serum (used as negative control) for free; ③. Rabbit polyclonal antibodies purified by Antigen Affinity;
Product Precision	1. CUSABIO guarantees the antibody purity > 95% confirmed by SDS-PAGE. 2. CUSABIO guarantees the antibody titer > 1: 64,000 confirmed by ELISA. 3. CUSABIO guarantees to offer Peptide-specific Dot blot. 4. CUSABIO guarantees 100% risk-free: No qualified antibody, no charge.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.