



P1 Antibody

Product Code	CSB-PA654400ZA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q39172
Immunogen	Recombinant Arabidopsis thaliana NADP-dependent alkenal double bond reductase P1 (P1) (1-345AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADPH-dependent oxidoreductase 2-alkenal reductase (AtAER) (EC 1.3.1.-) (EC 1.3.1.74) (NADP-dependent alkenal double bond reductase P1) (DBR1) (NADPH-azodicarbonyl/quinone reductase) (NADPH:2-alkenal/one alpha,beta-hydrogenase) (ALH) (P1-zeta-crystallin protein) (P1-ZCr), AER, P1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order
Gene Names	P1
Product Modify	Unmodified