




clfA Antibody

Product Code	CSB-PA753416ZA01SKW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GB45
Immunogen	Recombinant Staphylococcus aureus Clumping factor A (clfA), partial (228-558AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Clumping factor A (Fibrinogen receptor A) (Fibrinogen-binding protein A), clfA
Immunogen Species	Staphylococcus aureus (strain MSSA476)
Target Names	clfA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: clfA antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 52.1 kDa

Observed band size: 54 kDa

Product Modify	Unmodified
-----------------------	------------

Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.
--------------	---