



Sult1a1 Antibody

Product Code	CSB-PA022933ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52840
Immunogen	Recombinant Mus musculus Sult1a1 protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sulfotransferase 1A1 (ST1A1) (EC 2.8.2.1) (Aryl sulfotransferase) (Phenol sulfotransferase) (Phenol/aryl sulfotransferase) (mSTp1) (ST1A4) (Sulfokinase), Sult1a1, St1a1 Stp Stp1
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Sult1a1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.