



# Sult1a1 Antibody

<b>Product Code</b>	CSB-PA022933ZA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P52840
<b>Immunogen</b>	Recombinant Mus musculus Sult1a1 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mus musculus
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Sulfotransferase 1A1 (ST1A1) (EC 2.8.2.1) (Aryl sulfotransferase) (Phenol sulfotransferase) (Phenol/aryl sulfotransferase) (mSTp1) (ST1A4) (Sulfokinase), Sult1a1, St1a1 Stp Stp1
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	Sult1a1
<b>Product Modify</b>	Unmodified