



pyrG Antibody

Product Code	CSB-PA006176ZA01SMT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P65925
Immunogen	Recombinant Streptococcus pyogenes serotype M1 pyrG protein
Raised In	Rabbit
Species Reactivity	Streptococcus pyogenes serotype M1
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	CTP synthase (EC 6.3.4.2) (Cytidine 5'-triphosphate synthase) (Cytidine triphosphate synthetase) (CTP synthetase) (CTPS) (UTP--ammonia ligase), pyrG
Immunogen Species	Streptococcus pyogenes serotype M1
Lead Time	Made-to-order (12-14 weeks)
Target Names	pyrG
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.