



PB000185.00.0 Antibody

Product Code	CSB-PA307560ZA01EWN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7SI97
Immunogen	Recombinant Plasmodium berghei PB000185.00.0 protein
Raised In	Rabbit
Species Reactivity	Plasmodium berghei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-lactate dehydrogenase (EC 1.1.1.27)
Immunogen Species	Plasmodium berghei (strain Anka)
Lead Time	Made-to-order (12-14 weeks)
Target Names	PB000185.00.0
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.