



xylE Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA356892ZA01FFZ |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P06622 |
| Immunogen | Recombinant Pseudomonas putida xylE protein |
| Raised In | Rabbit |
| Species Reactivity | Pseudomonas putida |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Metapyrocatechase (MPC) (EC 1.13.11.2) (CatO2ase) (Catechol 2,3-dioxygenase), xylE |
| Immunogen Species | Pseudomonas putida (Arthrobacter siderocapsulatus) |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | xylE |
| Product Modify | Unmodified |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |