




ply Antibody

Product Code	CSB-PA312102ZA01FNF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04IN8
Immunogen	Recombinant Streptococcus pneumoniae serotype 2 ply protein (2-471aa)
Species Reactivity	Streptococcus pneumoniae serotype 2
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptococcus pneumoniae serotype 2 (strain D39 / NCTC 7466)
Target Names	ply

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: ply Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 66 kDa

Observed band size: 66 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.