



Ovomucoid Antibody

Product Code	CSB-PA356564ZA01CRB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05600
Immunogen	Recombinant Coturnix delegorguei Ovomuroid protein
Raised In	Rabbit
Species Reactivity	Coturnix delegorguei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ovomucoid (Fragment)
Immunogen Species	Coturnix delegorguei (Harlequin quail)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.