




# Trypsin/alpha-amylase inhibitor CMX2 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA674868ZA01TQN   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                         |
| <b>Uniprot No.</b>         | Q43691  |
| <b>Immunogen</b>           | Recombinant Triticum aestivum Trypsin/alpha-amylase inhibitor CMX2 protein (25-121aa) |
| <b>Species Reactivity</b>  | Triticum aestivum   |
| <b>Tested Applications</b> | ELISA, WB   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4      |
| <b>Purification Method</b> | Affinity-chromatography   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Product Type</b>        | Polyclonal Antibody   |
| <b>Immunogen Species</b>   | Triticum aestivum (Wheat)   |
| <b>Target Names</b>        | N/A   |

## Image



### Western Blot

Positive WB detected in: recombinant protein  
 All lanes: Trypsin/alpha-amylase inhibitor CMX2  
 Antibody at 1:1000  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/50000 dilution  
 Predicted band size: 33 kDa  
 Observed band size: 33 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.