



L3 Antibody

Product Code	CSB-PA361000ZA01HIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03253
Immunogen	Recombinant Human adenovirus C serotype 5 L3 protein
Raised In	Rabbit
Species Reactivity	Human adenovirus C serotype 5
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protease (EC 3.4.22.39) (Adenain) (Adenovirus protease) (AVP) (Adenovirus proteinase) (Endoprotease), L3
Immunogen Species	Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5)
Lead Time	Made-to-order (12-14 weeks)
Target Names	L3
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.