



BETVIII Antibody

Product Code	CSB-PA337344ZA01BSS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43187
Immunogen	Recombinant Betula pendula BETVIII protein
Raised In	Rabbit
Species Reactivity	Betula pendula
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding allergen Bet v 3 (Allergen Bet v III) (allergen Bet v 3), BETVIII
Immunogen Species	Betula pendula (European white birch?)(Betula verrucosa)
Lead Time	Made-to-order (12-14 weeks)
Target Names	BETVIII
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.